Granulox®

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PZN: 09505753

GB Instruction for use

Read all of this leaflet carefully because it contains important information for you. This medical device is available without prescription. However, you still need to use Granulox® carefully to get the best results from it.

- Keep this leaflet. You may need to read it again.
- Consult your doctor or pharmacist if you need more information or advice .
- You must contact a doctor if your symptoms worsen.
- If any of the side effects gets serious, or if you notice any side effects not listed in this leaflet, please tell your doctor or pharmacist.

In this leaflet:

1. What Granulox® is and what it is used for
2. What you need to know before you use Granulox®
3. What you need to know while using Granulox®
4. How to use Granulox®
5. Possible side effects
6. How to store Granulox®
7. Further information

1. What Granulox® is and what it is used for

Granulox® is a spray for use on chronic wounds. Granulox® is an innovative medical device for the treatment of chronic wounds, such as venous leg ulcer, arterial leg ulcer, mixed leg ulcer, diabetic foot ulcers, secondary healing of surgical wounds and pressure sores.

Chronic wounds are usually characterised by a below-average oxygen supply to the affected tissue. This has a major impact on wound healing. The precondition for use is appropriate medical treatment of the primary disease. In the case of chronic arterial wounds, the possibilities for improving blood flow must already have been exhausted.

Granulox® provides the wound with the required oxygen by means of diffusion. The active substance haemoglobin supplies the base of the wound externally with oxygen. The improved oxygen supply to the base of the wound supports wound healing.

2. What you need to know before you use Granulox®

Do not use Granulox®

- on infected wounds
- if you are pregnant.

No data are available to allow an assessment of these cases.

3. What you need to know while using Granulox®

- Granulox® should not be used simultaneously with locally effective medicines, such as antibiotics, as interactions have not yet been studied sufficiently.
- Disinfectants are known to impair the efficacy of Granulox®. Therefore, after application of a disinfectant, the area must be subsequently flushed thoroughly with a physiological solution.
- Thorough flushing with a physiological solution is necessary after use of proteolytic (enzymatic) debridements.

If your symptoms worsen or you experience any side effects, you must stop using the medical device and consult a doctor.

4. How to use Granulox®

Before you start to use Granulox®, a doctor or medically qualified person must have explained Granulox® treatment to you.

After sufficient debridement and cleaning of the wound, spray Granulox® onto the cleaned wound. Remove the protective cap from the spray container. Hold the nozzle about 5 to 10 cm from the wound and spray until the wound is completely coated with a thin film. 1 to 2 seconds are normally sufficient to cover a wound area of 2 x 3 cm. Repeat the application daily. After use, cover the wound with a breathable wound dressing. IMPORTANT: only use breathable wound dressings and bandaging materials.

After use, replace the protective cap on the spray container. It is not necessary to clean the nozzle as the product has been preserved. If the nozzle becomes blocked after several uses, it can be re-opened by cleaning it. Use sterile materials or disinfectants for cleaning. After cleaning, the first spray should be discarded, e.g. by spraying it onto cellulose.

Depending on the wound size, one pack of Granulox® is sufficient for a period of use of up to 30 days. After opening, Granulox® can be used until the expiry date.

5. Possible side effects

There are no known side effects of Granulox®. If your symptoms worsen or you develop side effects, please stop using Granulox® and consult a doctor.

6. How to store Granulox®

Granulox® should be stored in the refrigerator at temperatures from 2 °C to 8 °C. On the days it is being used, Granulox® can be stored at room temperature (max. 25 °C). If Granulox® is used on a daily basis (e.g. in clinical practice or home visits), the spray container can be stored continuously at room temperature until it is empty. The storage time at room temperature must not exceed a maximum period of 6 weeks. Store out of the sight and reach of children! Pressurised container. Store protected from sunlight and temperatures above 50 °C. Do not open forcibly or incinerate, even when empty. You may not use Granulox® after the expiry date given on the product and the outer packaging.

7. Further Information

What Granulox® contains

- The active ingredient is 10% carbonylated haemoglobin.
- The other ingredients are: 0.7% phenoxyethanol; 0.9% sodium chloride; 0.05% N-acetylcysteine; ad 100% water

What Granulox® looks like and contents of the pack

- Granulox® is a red, aqueous spray. One spray contains 12 ml.

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