## NANOGEN AKTIGEL INSTRUCTIONS FOR USE (English)

**Description:** Nanogen Aktigel is a bio-cellulose wound healing acceleration gel, which is derived from a unique plant extract. The product provides a variety of nutrients, enzymes and vitamins to the wound bed, the combination of which have been proven to accelerate wound healing. Nanogen Aktigel is proven to act at a cellular level within the wound. It requires the use of a suitable secondary dressing (GentilFoam) to prevent wound dehydration.

**Indications:** It can be used on most partial and full thickness wounds, including Leg Ulcers, Pressure ulcers, Diabetic ulcers, trauma wounds, complex surgical wounds and Epidermolysis Bullosa. (EB) It is particularly effective at kick starting chronic wounds which have become stuck in the inflammatory phase of healing. Suitable for application on all stages of wound healing, except wounds with more than 50% black necrotic tissue. Nanogen Aktigel is a 100% natural antimicrobial product and can be used on infected and critically colonized wounds.

**Contraindications:** None.

**Method of Usage:** Make sure that wound bed is cleaned prior to application with sterile saline. Spread the Nanogen Aktigel liberally across the base of the wound bed (to a depth 3-4 mm), but it must cover the whole wound, allow the gel to be absorbed into the wound bed (2-3 minutes). Once the Gel is placed on to the wound, cover with a suitable secondary dressing (GentilFoam) to prevent wound dehydration. It can also be used in conjunction with Nanogen Aktiv bio-cellulose matrix dressing to facilitate optimal performance of the Nanogen Aktiv product.

**Warnings:** 1. A suitable secondary dressing must be used. (GentilFoam) 2. The patient may experience a mild 'stinging' effect immediately after application; this is normal and will soon dissipate.

**Frequency of change:** Nanogen Aktigel should be reapplied each time the wound is re-dressed (the secondary dressing) until the wound is healed. The majority of the Nanogen Aktigel will be absorbed into the wound bed within a few minutes of application, but will continue to help stabilize the wound and create the correct environment for accelerated healing to occur.

**Storage:** Nanogen Aktigel dressings should be kept in a well-ventilated, dark environment with less than 80% relative humidity, between 10°C and 40°C temperature.

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