# OASIS<sup>®</sup>

## EXTRACELLULAR MATRIX

### Manufacture:

FF0079-01C  
Rev. 08/17  
Cook Biotech Inc.  
1425 Innovation Place  
West Lafayette, IN 47906, U.S.A.  
TEL: 765-497-3355  
FAX: 765-497-2361  
www.cookbiotech.com

### EC REP IC Representative:

Cook Ireland  
O'Halloran Road  
National Technological Park  
Limerick, IRLAND  
TEL: 353-413-3440  
FAX: 353-413-3444

### smith&hepner

Smith & HEPNER, Inc.  
TEL: 817-900-4000  
FAX: 817-900-4100

### ENGLISH

**INTENDED USE:**  
OASIS<sup>®</sup> EXTRACELLULAR MATRIX (ECM) is indicated for the management of wounds including:

- Partial and full-thickness wounds
- Pressure ulcers
- Venous ulcers
- Chronic vascular ulcers
- Diabetic ulcers
- Trauma wounds (abrasions, lacerations, skin tears)
- Second-degree burns
- Donor sites/grfts
- Draining wounds
- Surgical wounds (post-Mohs' surgery, post-laser surgery, podiatric, wound dehiscent)

OASIS ECM is supplied sterile in peel-open packages and is intended for one-time use. **CAUTION:** This product is intended for use by trained medical professionals.

**COMPOSITION:**  
OASIS ECM is derived using a process that retains the natural composition of matrix molecules such as collagen (Types I, III, IV, VI), glycosaminoglycans (hyaluronic acid, chondroitin sulfate A and B, heparin, and heparan sulfate), glycoproteins (fibronectin) and growth factors (FGF-2, TGF- $\beta$ ).

**CONTRAINDICATIONS:**  
This device is derived from a porcine source and should not be used in patients with known sensitivity to porcine material. This device is not indicated for use in third degree burns.

**PRECAUTIONS:**

- Do not re-sterilize. Discard all open and unused portions of the OASIS ECM.
- Device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- The device must be used prior to the expiration date.
- Discard device if mishandling has caused possible damage or contamination.
- OASIS ECM should not be applied until excessive exudate, bleeding, acute swelling, and infection is controlled.

**POTENTIAL COMPLICATIONS:** The following complications are possible. If any of these conditions occur, the device should be removed.

- Infection
- Chronic inflammation (Initial application of wound dressings may be associated with transient, mild, localized inflammation.)
- Allergic reaction
- Excessive redness, pain, swelling, or blistering

**STORAGE:** This device should be stored in a clean, dry location at room temperature.

**STERILIZATION:** This device has been sterilized with ethylene oxide.

**SUGGESTED INSTRUCTIONS FOR USING OASIS EXTRACELLULAR MATRIX**  
*These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.*

**NOTE: Always handle the OASIS ECM using aseptic technique.**

1. Prepare wound area using standard of care methods to ensure wound is free of debris and necrotic tissue. If necessary, surgically debride and treat the wound bed with standard of care techniques to ensure the wound shows granulation tissue and the edges contain viable tissue.
2. Cut the dry sheet to a size slightly larger than the outline of the wound area. If the wound is larger than a single sheet, then multiple sheets may be used. Overlap adjoining sheets to provide coverage of the entire wound.
3. Apply OASIS ECM in its dry state over the wound. Rehydrate the device with sterile saline or other isotonic solution. Alternatively, rehydrate OASIS ECM by placing it in a bowl of sterile saline or other isotonic solution for at least one (1) minute prior to use.
4. Place the edge of the device in contact with the intact tissue. Smooth the OASIS ECM into place to ensure the sheet is in contact with the underlying wound bed.
5. After application, use an appropriate, non-adherent, secondary dressing to maintain a moist wound environment. Change the secondary dressing as needed to maintain a moist, clean wound area. Frequency of secondary dressing change will be as needed.
6. Discard any unused portions according to institutional guidelines for biological waste. Reapply, as needed, if the OASIS ECM is no longer visible. Typically, reapplication is needed every 3-7 days until the wound is re-epithelialized.

**IMPORTANT:** As healing occurs, sections of the device may darken in color or form into a caramel-colored gel which may begin to peel from the area surrounding the wound. Any product that has lifted from the skin may be trimmed, taking care not to forcibly remove sections of the OASIS ECM that may adhere to the wound.

**EXTRACELLULAR MATRIX** This symbol means extracellular matrix.

**Fenestrated** This symbol means fenestrated

**Meshed** This symbol means meshed

**Bilaminat** This symbol means bilaminat

**Professional Sample** This symbol means professional sample

**REFERENCES:**

1. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med* 2007;18:537-543.
2. Hodde JP, Badyjak SF, Brightman AO, et al. "Glycosaminoglycan Content of Small Intestinal Submucosa: A Bioscaffold for Tissue Replacement." *Tissue Engineering*, 1999; 2: 209-217.
3. Data on File, Cook Biotech, Inc.

**REFERENCES: See English "REFERENCES"**  
© 2014 Cook Biotech Incorporated

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