SORELEX®**Productinformatie**

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SORELEX®

Product information

SORELEX® is a sterile medical device designed for covering and active healing of chronic wounds. It consists of an active layer and an absorbent core. Its comprehensive effect provides a favourable environment for moist wound healing. The antimicrobial effects of SORELEX® have a positive effect on the cleaning and decontamination of a wound during the initial healing phase, the active contact layer supports wound healing and the absorbent core absorbs the exudates from the wound.

Anti-adhesion properties and moist healing

Sodium hyaluronate is a linear, negatively-charged polysaccharide that is a natural part of the body (especially the extracellular space), therefore it is non-toxic and does not cause allergic reaction. Sodium hyaluronate has unique physical-chemical properties, a high hydrating ability and excellent anti-adhesive properties. Thus, a layer of sodium hyaluronate on the contact layer of SORELEX® provides an appropriate environment for the wound's natural healing process.

Antimicrobial efficacy

Octenidine - N, N' - (1, 10-decanediyl-di-1 [4 h]-pyridinyl-4-ylidene) bis - (1-octanamine) dihydrochloride – is an antiseptic used to prevent infections that slow down the wound healing process. Octenidine has a broad spectrum of antimicrobial activities against many gram-positive and gram-negative bacteria, fungi and plaque-producing organisms. Its onset is faster than other antiseptics (e.g. silver compounds, chlorhexidine, PVP-iodine, triclosan and polyhexanide) and it has a wider range of effects even at very low concentrations. Octenidine is not absorbed by intact skin or wounds – thus it does not have an effect on the entire body making it toxicologically harmless even during repeated use.

Ingredients

Sodium hyaluronate, octenidine dihydrochloride

Indications

SORELEX® is used to treat wide spectrum of non-healing chronic wound. SORELEX® may be used on both infected and stalled wounds. Primary indications are leg ulcers, diabetic foot ulcers, pressure ulcers and dehisced surgical wounds.

Application

Observe the principles of good hygiene and maximum cleanliness when applying the product.

The contact layer of SORELEX® is applied directly to the wound, the coloured backing layer is aligned away from the wound. When applied to a wound with very little secretion it is necessary to moisten the contact layer with a small amount of sterile saline solution or water (about 10 ml). Moistening is not necessary for other wounds. Depending on the degree of exudation (e.g. absorption protection) a suitable secondary dressing should be chosen to apply over SORELEX® or retain with a bandage fixed to the wound. SORELEX® may be used up to 5 days, however it is recommended to change SORELEX® at 2-3 day intervals depending on the wound's condition. For highly exuding wounds the SORELEX® wound contact layer may be left undisturbed on the wound bed and the secondary dressing changed.

During dressing changes it is necessary to pay attention to the healing edges of the wound, which are covered by new, soft skin. It is possible to apply an appropriate salve to prevent the SORELEX® clinging to the new epithelium (skin) on the edges of the wound. If SORELEX® dries out, to prevent adherence to new formed tissue, soaking the cover will facilitate its removal from the wound when changing. If the SORELEX® does not significantly overlap the edges of the wound and is fully soaked with secretion, the cover's textile layers may be disconnected allowing a more careful removal from the wound (first remove the top absorbent part, then carefully remove the light contact textile).

Attention

If there is any deterioration in the wound's state or an unexpected condition, contact your doctor.

Contraindications

Currently there are no known contraindications for the local use of hyaluronate sodium. The octenidine dihydrochloride contained in the product may invoke an unfavourable reaction in individuals with hypersensitivity. Therefore SORELEX® is contraindicated for patients who are hypersensitive to octenidine.

Side effects

Sodium hyaluronate is well tolerated. Currently there are no known side effects after local administration.

A sensation of burning and heat at the point of application were occasionally observed. Octenidine dihydrochloride may invoke an unfavourable reaction in individuals with hypersensitivity.

Any adverse reactions or other unusual reactions should be discussed with your doctor.

Any serious adverse reactions to SORELEX® should be reported to the manufacturers:

Contipro a.s., Dolní Dobrouč 401, 561 02 Dolní Dobrouč

tel: +420 800 331 157, +420 465 520 035, fax: +420 465 524 098, e-mail: registration@contipro.com

Interactions

The product should not be combined with preparations containing anionic surfactants (e.g. detergents and cleaning products), with which it may form insoluble complexes. This reduces its effectiveness. The efficiency of SORELEX® is also reduced in combination with products containing iodine, PVP-iodine or other oxidising substances.

Warning

Only for external use

Do not use a product with damaged packaging.

Do not use after the expiry date indicated on the package.

Do not use if the contact layer is damp (this may be caused by gelation or adhesion of the layer to the package).

Keep out of reach of children.

Package contents

1/ 5/ 10 pcs SORELEX®, information leaflet.

Shelf life

24 months

Use immediately after opening the primary container

Storage

Store in a dry place at a temperature of 2 °C to 25 °C. Protect from light and frost.

Manufacturer

Contipro a.s.

Dolní Dobrouč 401

561 02 Dolní Dobrouč

Česká republika

Symbols used

Consult instructions for use

Temperature limit

Batch code

Use by date

Sterilized using ethylene oxide

CE mark

Do not use if packaging is damaged

Manufacturer

Date of manufacture

Keep dry

Single use

Do not re-sterilize