BioMonde

Wat is

# LDT Madentherapie?

Making healing possible

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1. Cazander G et al. 2013. Multiple actions of *Lucilia sericata* larvae in hard-to-heal wounds. *Bioessays*; 35: 0000–0000.
2. Markevich Y, McLeod-Roberts J, Mousley M, Melloy E. 2000. Maggot therapy for diabetic neuropathic foot wounds: a randomized study. *Diabetologia*; 43: A15. Suppl.
3. Mudge E, Price P, Walkley N, and Harding K. 2014. A randomized controlled trial of larval therapy for the debridement of leg ulcers: results of a multicenter, randomized, controlled, open, observer blind, parallel group study. *Wound Repair & Regeneration*; Jan-Feb, 22(1); 43-51.
4. Cowan LJ, Stechmiller JK, Phillips P, Yang Q, et al., 2013. Chronic wounds, biofilm and use of medicinal larvae. *Ulcers*; 487024.
5. Iková H et al. 2013. Growth and Survival of Bagged *Lucilia sericata* Maggots in Wounds of Patients Undergoing Maggot Debridement Therapy. *eCAM*; 192149.
6. Blake, F. 2007. The biosurgical wound debridement: Experimental investigation of efficiency and practicability. *Wound Rep Reg*; 15(3): 756-761.
7. Opletalová, K. et al.: Maggot Therapy for Wound Debridement. A Randomized Multicentre Trial. *Arch. Dermatol.*, published online December 19, 2011
8. Gilead L, Mumcuoglu KY, Ingber A. 2012. The use of maggot debridement therapy in the treatment of chronic wounds in hospitalized and ambulatory patients. *J Wound Care*; 21(2):78–85.
9. Bültemann. 2015. Anleitung zur Applikation medizinischer Fliegenlarven (*Lucilia sericata*) bei der Behandlung chronischer Wunden. *Wund Management* 2015; 9 (3), S. 100-106
10. Schultz G. et al.: Wound bed preparation: a systematic approach to wound management. *Wound Repair and Regeneration* 11(2), 2003, Supplement
11. Horobin AJ, Shakesheff KM, Pritchard DI (2005) Maggots and wound healing: an investigation of the effects of secretions from *Lucilia sericata* larvae upon the migration of human dermal fibroblasts over a fibronectin coated surface. *Wound Repair Regen* 13(4): 422–33.
12. Bexfield A, Bond AE, Morgan C et al (2010) Amino acid derivatives from *Lucilia sericata* excretions/secretions may contribute to the beneficial effects of maggots therapy via increased angiogenesis. *Br J Dermatol* 162(3): 554–62. doi: 10.1111/j.1365-2133.2009.09530.x
13. van der Plas MJA, Baldry M, Van Dissel JT, et al (2009) Maggot secretions suppress pro-inflammatory responses of human monocytes through elevation of cyclic AMP. *Diabetologia* 52: 1962-70.
14. Bennett H. et al. 2013. Cost-Effectiveness of Interventions for Chronic Wound Debridement: an Evaluation in Search of Data. *Wounds UK*; 9(4): 3-11. Suppl.
15. Mumcuoglu, K et al. 2001. Destruction of Bacteria in the Digestive Tract of the Maggot of *Lucilia sericata* (Diptera: Calliphoridae). *Journal of Medical Entomology*, Volume 38, Issue 2, pp. 161–166,
16. Cazander G., et al. 2009. Maggot therapy for wound healing: clinical relevance, mechanisms of action and future prospects. *Journal Wound Technology*, 5: pp.18-23
17. Cazander G., et al. (2010) Maggot excretions inhibit biofilm formation on biomaterials. *Clin Orthop Relat Res*, 468: pp. 2789-2796.
18. Nigam, Y. 2013. LDT: An economic, scientific and clinical evaluation; Evidence for larval debridement therapy in wound cleansing and healing. *Wounds UK*; 9(4) Suppl.
19. Tian X et al. 2013. Maggot debridement therapy for the treatment of diabetic foot ulcers: a metaanalysis. *Journal of Wound care*; 22(9): 462-9.
20. Steenvoorde P, Jacobi CE, Van Doorn L, Oskam J. 2007. Maggot debridement therapy of infected ulcers: patient and wound factors influencing outcome – a study on 101 patients with 117 wounds. *Ann R Coll Surg Engl*; 89(6): 596–602.
21. Campbell N. 2014. A Retrospective, Quality Improvement Review of Maggot Debridement. *Therapy Outcomes in a Foot and Leg Ulcer Clinic. OWM*; 60(7): 16-25.
22. Armstrong et al.. 2005. Maggot therapy in "Lower-Extremity Hospice" Wound care. *Journal of the American Podiatric Medical Association* 95 (3), 254 - 257.
23. Sherman R. 2002. Maggot versus conservative debridement therapy for the treatment of pressure ulcers. *Wound Repair Regen*; 10; 208-14.

**BioMonde Ltd.**  
Units 2–4 Dunraven Business Park  
Coychurch Road  
Bridgend  
CF31 3BG  
United Kingdom

Telephone: +44 (0) 345 230 1810  
Fax: +44 (0) 1656 668 047  
Email: [enquiries@biomonde.com](mailto:enquiries@biomonde.com)

[www.biomonde.com](http://www.biomonde.com)

Zoobiotic Limited, trading in the UK as BioMonde.