

Cover

skin protectant

For the management and prevention of IAD, intertrigo and other moisture lesions



Cover is a topical dimethicone-based barrier ointment

- ✓ Suitable and safe for both intact and damaged skin
- ✓ Reduces the risk of skin damage due to friction and shear forces
- ✓ Helps manage skin irritations like redness and chapping associated with pressure ulcers (up to stage 2)
- ✓ Can be used on a daily basis or as often as required
- ✓ Transparent, for better assessment of skin conditions
- ✓ Good adhesion to both moist and dry skin, yet easy to apply and remove
- ✓ Zinc oxide free, therefore also suitable for patients with radiotherapy
- ✓ Fragrance-free
- ✓ Vegan



Cover

skin protectant



Cover bridges the gap between continence and wound care

Incontinence is a major problem for healthcare systems worldwide. With life expectancies in the Western world still on the increase, the problem will only get worse. Roughly half of all institutionalized elderly patients suffer from incontinence; urinary, faecal or both.

The results are all too familiar to care professionals everywhere: skin irritation, reddening, inflammation, infections and excoriations, all symptoms of incontinence associated dermatitis (IAD).

The strong association between IAD and pressure ulcers underscores the need for integrated incontinence and pressure ulcer prevention in every care setting.

Cover has been specifically designed for the daily skin care of incontinent patients. Therefore it represents a welcome contribution to the cost-effective prevention and treatment of moisture lesions, IAD and intertrigo.

Cover fulfils three important roles:

1. Cover helps to prevent the development of IAD, moisture lesions, pressure ulcers and intertrigo.
2. Cover helps to protect the healthy, irritated and/or broken skin by creating a barrier.
3. Cover helps to promote recovery by providing protection and an optimum moist environment.

Why Cover is cost-effective

Cover offers a practical and cost-conscious solution for the prevention and care of incontinence associated dermatitis (IAD) and other forms of moisture associated skin damage (MASD), including intertrigo and early-stage pressure ulcers. Cover supports both preventive skin care and the management of superficial skin damage.

- ✓ Prevents and contains early-stage IAD and MASD, reducing the need for more intensive interventions.
- ✓ Minimizes nursing time by enabling fast, uncomplicated daily application.
- ✓ Eliminates the need for secondary wound dressings, thanks to its protective barrier.
- ✓ Two practical tube sizes support tailored use.

Indications

Cover is indicated for the care of moisture-associated damage of intact or injured skin as a result of frequent incontinence, salivation, perspiration and exudation. It is also indicated for preventing and managing skin irritations like redness and chapping associated with pressure ulcers (up to category 2) and intertrigo.

Contra-indications

Do not use the product when the patient is known to be sensitive to one of the ingredients of the ointment (see Ingredients). Hypersensitivity to one of the ingredients of the ointment can cause redness, irritation and possibly light pain.

Instructions for use

First thoroughly cleanse the entire skin area with a suitable cleanser. Wash or cleanse hands, remove ointment from tube and apply a liberal, thick coating of Cover. The treatment may be repeated as frequently as is necessary. Always read the package insert of this product.

Packaging

Cover is available in tubes of:
60 gr (11040)
120 gr (11050)

CE marking

Cover is a Class I medical device and manufactured by: Bfactory Health Products B.V., The Netherlands.

Ingredients

Dimethicone, Glycerin, Propylene Glycol, Polyethylene Glycol 400, Polyethylene Glycol 4000, Povidone, Cetostearyl Alcohol.

Distributor

BAP Medical B.V.
The Netherlands
info@bap-medical.com
© 2026 BAP Medical B.V.

